

CLIA Waived

www.hemosense.com

1.877.436.6444 x5

INRATIO®

Why more Healthcare Professionals
are Choosing the...

INRatio® PT/INR Monitoring System

For consistent PT/INR results in under 2 minutes using
one drop of fresh whole blood from a fingerstick.

QC
you can
see!



2 Level On Board
Quality Controls
report **quantitative**
results you can see

Normal and therapeutic quantitative controls
that test integrity across the full reporting range
are performed side by side with the PT test, using
the same drop of blood, ensuring the reliability
of each and every test. No additional costly, time
consuming external control tests or materials
are necessary.



Individually Wrapped
Test Strips Do Not
Require Refrigeration

Store at room temperature for up to 12 months.
No need to wait for strips to warm up.
Save precious refrigerator space.



Small Sample Size

One unmeasured drop (≈15 µL) of fresh whole
blood from a fingerstick applied directly to the
strip performs both the test and quality controls.
No special blood collection devices required.



inverness medical

INRatio® is a registered trademark of the Inverness Medical group
of companies. © 2007 Inverness Medical. All rights reserved.

INRATIO

1.8.08 12:48PM

INR=2.6

PT=25.8

QC1=10.0

QC2=30.0

STRIP CODE=KC9ZC

menu

OK

^

v

INRATIO®

Actual Size

circle #108

Send patients home instead of to the lab!
Fast test results allow immediate consultation
for warfarin dosage changes. Efficiency for
the practice, convenience for patients!